





CLOUD

BEDHEAD UNIT BH-MEX005 for ICU



Spesification :

- Model : Bed head unit Full Oval for ICU
- Dimension : 1500mm x 215mm x 60mm.
- Alluminium Thickness Construction : 1.8mm
- Material : Alluminium alloy.
- Weight : 10.5 Kg

Accessories :

- 8x Electrical Socket.
- 1x on/off switch Light.
- 1x Reading LED Light.



 SUKET FR.03.02/VA/6366/2020



PT. AMSAR PRIMA MANDIRI

Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340
www.amsarmedical.com info@amsarmedical.com

Telp. (021) 75674728