

**KL-T-IB**

**Electrical Single Arm Surgical Pendant KL-T-IB** is a new and practical equipment used in modern clean operation rooms and general operation rooms.

## Spesification

- Single Arm : 1000
- Horizontal Rotation
- Electric moving up and down
- 2 tray
- 2 IV Pole
- 2 Baskets (European Standard)
- Gas Outlet (O, O, A4B, A7B, V)
- High Carrying Capacity
- Electric Driven
- Pneumatic Lifting

## Basic Configuration

- one discharge device (German style) for compressed air, oxygen, nitrous oxide, suction gas and anesthetic gas respectively, four power outlets, one network outlets. (Optional gases)

## Features

- The control cabinet of the tower crane adopt the imported high-strength aluminum alloy of patent technology for the one-off extrusion molding, and the surface undergoes the primary oxidation treatment.
- Aluminum molding tray with non-slip texture.
- Mechanical friction damping brake.
- The imported electrical machine is adopted, and the rising and rotation function of the electric horizontal cantilever can accurately position without excursion.
- The gas pipeline, power supply and computer communication line are separately arranged without interference.
- The imported German standard gas terminal (over 20,000 of inserting and pulling out) is adopted.

**SUKET FR.03.02/VA/6131/2020**