

Keling KL-LED · STZ5 SERIES

Shadowless Operating Lamps



Jiangsu Keling Medical Appliances Co., Ltd.
Shanian Industrial zone, Gaoyou City, Jiangsu Province



KL-LED.STZ5 series LED shadowless operating lamp is suitable for various operation needs. It is a modern ideal surgical lighting equipment operating room and purify the operating room.

Features

- Adopt imported LED cold light source, white the ultra-thin optical lens.
- Adopt digital means for stepless regulation of LED brightness.
- Adjust the diameter of light spot and color temperature.
- Removable handle sheath, easy to clean and disinfect.
- Flower-shaped super thin design light head.
- With Endoscopy model.
- Titanium alloy balance arm, no welding, never rust.
- Nano coating of heat sink has created an excellent heat dissipation effect.

Specification

MODE	KL-LED-STZ5
Diameter of lighthouse (cm)	Combined type
Illumination at 1m (lux)	80000-160000
Light field diameter at 1m (cm)	10-28
D10 (mm)	180±40
Diameter of spot D50 (mm)	110±35
L1 + L2 Depth of illumination (cm)	150
Color rendering index (Ra)	100≥Ra≥85
Optional color temperature (K)	6700≥Tc≥3000(K)
Number of LED bulbs	90/90
Average bulb life (hours)	50000
Heat-to-light ratio (mW/m ² lx)	< 6
Working distance (cm)	70-140
Power Supply	100-240V 50/60HZ
Input power (VA)	180
Ceiling height for installation (cm)	270-320

KEMENKES RI AKL 21603025921



Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340

PT. AMSAR PRIMA MANDIRI

www.amsarmedical.com info@amsarmedical.com

Telp. (021) 75674728