

# SOLIDA NITROUS OXIDE OUTLET CIG TYPE



## Spesifikasi

Outlet CIG Type

Available: Nitrous Oxide(Blue)

Material: Zinc, stainless steel, and plastic

Dimension: Width 5 Inch (127 mm) and height, 5.5 inch (140MM)

 **SUKET FR.03.02/VA/05624/2021**

**PT. AMSAR PRIMA MANDIRI**

Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340

[www.amsarmedical.com](http://www.amsarmedical.com) [info@amsarmedical.com](mailto:info@amsarmedical.com)

Telp. (021) 75674728