


Solida

SOLIDA 02 FLOWMETER - BS F25, 0-5 LPM



Features ↘

- ▷ Product Code : 09-1-02005.
- ▷ 0-5 LPM flow rate, the maximum flow ranges limited to 60 LPM (For Pediatric).
- ▷ Back inlet connection.
- ▷ Body made of chromed-brass.
- ▷ Flow control knob is in front of the item.
- ▷ Level indicator is green and flow indicator ball is stainless steel.
- ▷ CE marked.
- ▷ Pop-off valve on cover automatically open at 2 psi as the outlet line occlude.
- ▷ Whistle alerts the user to deal with it immediately.
- ▷ Recessed connector with screw cap protects the connector from accident breakage.
- ▷ 360° DISS swivel connection for easy installation.
- ▷ DISS female connectors in metal construction for maximum strength and durability.
- ▷ Polycarbonate shatter-resistant bottle for maximum durability and service life.
- ▷ Maximum and minimum water level markings clearly defined.
- ▷ High performance diffuser provides small bubbles for quiet operation with high humidify output
- ▷ Produce very low noise helps the patient rest.
- ▷ Back inlet connection.
- ▷ Body made of chromed-brass.
- ▷ Flow control knob is in front of the item.
- ▷ Level indicator is green and flow indicator ball is stainless steel.



Spesifikasi ↘

- ▷ The listed Flowmeter are accuracy tested at 10% tolerance.
- ▷ The listed Flowmeter Are calibrated at specified inlet pressure, 21°C standar atmospheric pressure.
- ▷ Calibration and color may be varied by customer requirment.
- ▷ Pressure other than indicated may affect flometer accuracy .
- ▷ Gas use appropriate gas adapter/ connector for connecting to the proper gas source .
- ▷ Teperatur transport and storage -20°C - 60°C.

Solida



 KEMENKES RI AKL 20403918270



PT. AMSAR PRIMA MANDIRI
Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340

www.amsarmedical.com info@amsarmedical.com

Telp. (021) 75674728