

Solida

TWIN O2 FLOWMETER



TWIN O2 FLOWMETER with O2 ADAPTER OHMEDA TYPE

Spesification

- ▷ 0-15 LPM flow rate, flood flow ranges from 60 LPM, calibrated at 21° C.
- ▷ Flowmeter body made of chromed-brass.
- ▷ Pressure control knob is in the front side.
- ▷ Flow indicator ball is stainless steel.
- ▷ USA color code.
- ▷ Flowmeter is CE marked.
- ▷ Attached with Manufacturer Calibrated/-tested.
- ▷ Y connector made of chromed-brass.
- ▷ Y connector is 1/8" NPT F, customizable with different adapter styles. (Ohmeda style, Silberman style, JIS style, Chemetron style, Puritan style, Drager style, MGI style, Etc.)
- ▷ Made in Thailand.
- ▷ ISO 13485: 2016.

Solida



 KEMENKES RI AKL 20402918270



PT. AMSAR PRIMA MANDIRI

Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340

www.amsarmedical.com info@amsarmedical.com

Telp. (021) 75674728