




Toilet Grab Bar Series



HS-018B-3



Specification:

Material: Nylon an metal composite

Configuration: Ordinary fittings

Colour: white, yellow

Size: 600*135 mm

FEATURES :

Unique anti-bacterial technology to improve the health enviroment.

Floating point and anti-slip design make it more safe and comfortable to grip.

Self-extinguishing material has high melting point, which is not combustion-supporting.

Anti-static and non-dust absorption, easy to clean, wear-resistant, water-resistant, acid and alkali resistant.

Easy to install, flexible combination, suitable for different installation condition.



Suket: FR.03.02/VA/04034/2021



PT. AMSAR PRIMA MANDIRI

Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340

www.amsarmedical.com - info@amsarmedical.com

Telp (021) 75674728